# Noviomagus Revision Mesh

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## EN Instructions for use

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## NO Bruksanvisning

## SV Bruksanvisning

# Noviomagus Revision Mesh

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# Noviomagus Revision Mesh

## EN INSTRUCTIONS FOR USE

Read these instructions before use

### PRODUCT INFORMATION

The Noviomagus Revision Mesh is a stainless steel mesh for orthopaedic use. The Noviomagus Revision Mesh is available in the types shown on the front page: Universal Flat (1), Acetabular Rim Large (2), Acetabular Rim Small (3), Proximal Femur (4), Acetabular Medial Wall Small (5) and Acetabular Medial Wall Large (6) and is designed for the reconstruction of the acetabulum and femur in combination with Bone Impaction Grafting (BIG) during a hip reconstruction. For different indications specific types of meshes are designed:

#	Type of Mesh	Indication
1	Universal Flat	General defects of the femur and the acetabulum
2	Acetabular Rim Large	Defects of the rim of the acetabulum
3	Acetabular Rim Small	Defects of the rim of the acetabulum
4	Proximal Femur	Defects of the proximal part of the femur
5	Acetabular Medial Wall Small	Defects of the medial wall of the acetabulum
6	Acetabular Medial Wall Large	Defects of the medial wall of the acetabulum

### MATERIAL

Implant grade austenitic stainless steel.

### INDICATIONS

The Noviomagus Revision Mesh is solely meant to contain bone graft during BIG techniques to reconstruct bone. The mesh is not for load bearing use.

### KONTRAINDICATIONS

The Noviomagus Revision Mesh must not be used when the patient has:

- absolute contraindications:
  - infection,
  - osteomyelitis,
  - serious osteoporosis at the place of fixation of the mesh,
  - malignancy in the surgical area,
  - neurological or mental disorder leading to unwillingness and/or disability to follow instructions.
- relative contraindications:
  - metabolic disorders which may impair bone function,
  - osteomalacia,
  - distant foci of infection which may cause haematogenous spread to the implant site,
  - confirmed or suspected metal sensitivity.

### INSTRUCTIONS FOR USE

- Determine the need for a revision mesh to contain bone graft during a hip reconstruction using BIG. Determine the type of mesh.
- Open the packaging.
- Determine whether the anatomic shaped Noviomagus Revision Mesh or special mesh screws fits the patient's anatomy. If necessary, cut the Noviomagus Revision Mesh to the right size with a steel-scissor. Take into account the sharp edges caused by cutting. If necessary, bend the mesh by hand into the desired shape.
- Place the Noviomagus Revision Mesh on position. Fixate the Acetabular Rim and Acetabular Medial Wall meshes with bone or special mesh screws with a minimal head diameter of 4 mm. Fixate the Proximal Femur mesh with cerclage wire. Special elevatable loops are integrated in the Acetabular Rim and Universal Flat mesh. The Universal Flat mesh can also be fixated with cerclage wire as well as screws.

### CAUTION

- The Noviomagus Revision Mesh must not be used in a microbiological contaminated site.
- This product is only to be handled and/or implanted

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All rights reserved. Printed in the Netherlands.