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### APPLICATION

This instrumentation consists of reusable manual instruments intended for use during surgical procedures.

### UTILIZATION

- Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation.
- If applicable, trial components should be used for size determination, preparation evaluation, trial reduction and range of motion evaluation, thus preserving the integrity of the actual implants and their sterile packaging.
- If applicable, available radiographic templates should be used to assist in the preoperative prediction of component size and type.
- A Surgical Technique can provide additional procedural information.

### WARNINGS

- Care should be taken not to cut through surgical gloves when handling any sharp-edged device.
- Due to differing parameters, varying tolerances, different materials and manufacturing specifications used by different manufacturers, this instrumentation should not be used to implant any other manufacturer's components. Any such use will negate the responsibility of Spierings Orthopaedics BV for the performance of the resulting implant.

### ADVERSE EFFECTS

- Malfunctions of instruments including, but not limited to, weld fracture, thread gauling and material fatigue, can occur in a small percentage of cases.
- Serious complications may be associated with any surgical procedure. These complications include, but are not limited to: genitourinary disorders;

gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders including emboli; myocardial infarction or death.

- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.
- Intraoperative fissure, fracture, or perforation of the bone can occur due to numerous factors including the presence of defects, poor bone stock, the use of the surgical instruments and the impaction of the components into the preparation.

### STERILIZATION

- The instrument is **not** sterile.
- Remove all instruments from their packaging prior to sterilization.
- All parts must be cleaned and sterilized before use.
- Spierings Orthopaedics BV Instruments are recommended to be cleaned and sterilized by the steam autoclaving procedure regularly used in the hospital. A suggested method is described in the document 'Reprocessing of reusable surgical instrumentation of Spierings Orthopaedics BV according to ISO 17664'.

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Spierings Orthopaedics B.V.