

Surgical Technique of the Bone Harvesting System

Application: This device is solely meant for harvesting bone chips (autograft) by means of drilling.




ASSEMBLY

Assemble the Bone Harvesting System according to the instructions below and the instruction figure (next page):

1. Choose the appropriate size Bone Harvesting Drill.
2. Place a *Bone Harvesting Drill* in the *Ratchet T-Handle* using the AO-coupling or a low RPM (up to 200 RPM) power drill using the upper part of the AO- coupling.
3. Put the ratchet T-Handle in position "R"

SURGICAL TECHNIQUE

1. After assembly, drill to a predefined approximate depth by rotating the drill clockwise.
 - Depth markings on the *Bone Harvesting Drill* have increments of 10 mm.
 - Size of the grafts can be varied from fine to coarse by applying more pressure to the drill while drilling and/or by decreasing the rotational speed of the drill.
2. Remove the *Bone Harvesting Drill* from the *T-Handle* or power drill.
3. Hold the rear end opening of the *Bone Harvesting Drill* above a bowl.
4. Push the harvested bone grafts through the *Bone Harvesting Drill* using the *Push Rod* (see instruction figure next page). Slight tapping on the handle of the *Push Rod* may be required to loosen bone grafts in the *Bone Harvesting Drill*.
5. Place the instruments back in the *Sterilization Tray* for cleaning and sterilization.




Warnings

- Be careful for the sharp cutting edge of the *Bone Harvesting Drill*.
- Be careful not to damage the cutting edge of the *Bone Harvesting Drill* when using the *Push Rod* to remove the bone grafts.

Remarks

- For instructions for use see document "Instructions for use of reusable surgical instrumentation" (IFU-MOI001).
- For reprocessing of reusable medical devices see document "Reprocessing of reusable medical devices" (REP-MOI005).

INSTRUCTION FIGURE



STERILIZATION TRAY

